



Burnshield

Burn Trauma First Aid Dressing

Life Saving Dressing for immediate use on 1st, 2nd & 3rd degree Burns and Scalds

PRINCIPAL MANUFACTURERS : LEVTRADE INTERNATIONAL PTY LTD ,
SOUTH AFRICA



Essential First Aid Dressing for :

- Fire Fighters , Rescue Teams
- Ambulance, Paramedic
- Emergency medical Services
- Hospitals, Clinics
- Police , Crowd Control squad
- Hotels, Restaurants
- Institutions, Laboratories
- Industrial Establishments
where burns are a hazard

NOW AVAILABLE

CLINICAL TEST EFFICACY REPORT FROM “SAFDARJUNG BURNS HOSPITAL, DELHI”
Asia's Largest Burns Hospital

SALIENT FEATURES:

- Provides immediate physical protection to the burn victim. Can be applied on all types of burns.
- Gives instant pain relief to a victim thru its cooling , moistening & trauma reducing properties.
- Soaked in non-irritant sterilised gel to eliminate the risk of infection from airborne bacteria.
- Instantly stops burn progression.
- Dressing is non adhesive and can be easily removed without causing pain to the victim.

DESCRIPTION

Sterile , individually packed Dressing made of special, soft elastic poly foam sheet for retaining - holding water based medicinal Gel . The Hydro Gel is non-toxic, non irritant, non staining, translucent off white colour in appearance & is composed of 96 % water with balance Melaleuca alt, emulsifiers etc

Available in various Sizes : 10 cm x 10 cm , 20 cm x 20 cm , 20 x 45 cm , 60 cm x 40 cm &
Easy Care Utility Packs too

Certifications / Reports

- Directorate General of Health Services, Ministry of Health , India.
- Safdarjung Burns Hospital , N.C.T. New Delhi
- Registrar of Medicines, South Africa for Distribution in India
- ISO 9001 & 13845 Quality Assurance Certification for OEM
- CE certificate for the OEM
- Sterility Test Certificate for every batch / Lot of Dressing