

BURNSHIELD®

BURNSHIELD®, when applied to burn wounds, preferably as soon as possible after the burn, has the tremendous ability to absorb and dissipate heat, providing relief, minimising shock and skin damage, halting the burn process, and physically protects against further contamination. BURNSHIELD®'s cooling and soothing effect results from its specially formulated Hydrogel which is non-toxic and non-irritant.

The active ingredient is Melaleuca alternifolia. The sterile dressing keeps the wound clean while

the loss of body fluids is kept to a minimum. BURNSHIELD® facilitates removal of burnt clothing. Easy to apply, non-adhesive and easy to remove without further trauma to the patient, allowing for quick medical diagnosis. Safe for use on children.

BURNSHIELD® can be used as an emergency dressing for superficial to full thickness burns. Any burn, whether received from open flames, scalding from hot water, steam, gas, chemicals or electrical appliances can be treated immediately with BURNSHIELD®.

For immediate use on burns and scalds. BURNSHIELD® blankets and dressing features are:

- Sterile - bacteria static
- Suitable for use on most burns and scalds
- Non-toxic
- Non-irritant
- Physical protection against contamination
- Cool and moistens the burnt area
- Minimises skin damage
- Inhibits the burn's progression
- Non-adherent and easy to remove
- Covered in water soluble gel
- Elastic and able to stretch and flex
- Versatile and innovative
- Portable
- Available in a variety of sizes
- An effective after-sun care in the form of Solace
- Stable gel structure



BU/095



BU/092




BU/093



BU/094



BU/091



BU/099



BU/086

BU/098



BU/097

BURNSHIELD TRAUMA BURN DRESSINGS

- Sterile trauma burn dressing
- Simple to apply and remove without any discomfort
- Open cell foam dressings, individually packaged and easy to open
- Contains hydrogel, minimises trauma
- The packet for BU/092 can also be used as an emergency burn bag - used in conjunction with the dressing
- Composed of 96% water and 4% tea tree oil and emulsifiers
- Simple to use - no special training required
- Provides instant cooling and pain relief
- Provides physical protection from infection
- Safe for use on children
- The dressing is non adherent and non toxic
- Ideal for pre-hospital burns care

- BU/091 Burnshield 100 x 100mm - single |
- BU/092 Burnshield 200 x 200mm - single |
- BU/093 Burnshield 600 x 400mm - single |
- BU/094 Burnshield 25 x 500mm - digit |
- BU/095 Burnshield 50 x 1000mm - limb |
- BU/098 Burnshield hydrogel 125ml |
- BU/086 Burnshield hydrogel 50ml |
- BU/099 Burnshield hydrogel burn blott - single |
- BU/097 Burnshield Face Mask 200 x 450mm |



EASycARE BURN KIT

- Suitable for small/medium burns
- Ideal for home, clubs, & DIY
- Supplied in a blue plastic box
- Contents:
 - 1 x Blue Plastic Box
 - 2 x 10 x 10cm Burnshield Dressing
 - 1 x 20 x 20cm Burnshield Dressing
 - 3 x Burn Gel in Sachet 1/8 oz.
 - 3 x Adhesive Plasters
 - 1 x 2.5cm x 5m Roll of Tape
 - 1 x Latex Gloves Pair
 - 1 x Scissors
 - 2 x 5cm Conforming Bandage

BU/087 Burnshield EasyCare Burn Kit



BU/087

BURNSHIELD CONTOUR KITS

- Suitable for large area burns, especially on the torso & back
- Each dressing is saturated with a sterile water based gel which moisturizes the burn area, cools and soothes the skin thereby minimises trauma
- Ideal for EMS and hospital
- Contour dressings available in 1m x 2m & 1m x 1m sizes
- Kits available in 2 sizes:

BU/073 Contains:

- 2x 1m x 1m Burnshield Contour Dressing
- 2x Hydrogel Bottle - 125ml
- 1x Soft Dressing 60mm x 5m

BU/074 Contains:

- 1x 2m x 1m Burnshield Contour Dressing
- 2x Hydrogel Bottle -125ml
- 1x Soft Dressing 60mm x 5m

BU/071 Burnshield Contour Dressing 1m x 1m

BU/072 Burnshield Contour Dressing 2m x 1m

BU/073 Burnshield Contour Kit incl. 2 x 1m x 1m

BU/074 Burnshield Contour Kit incl. 1 x 2m x 1m

BU/073



BU/071



BU/072

BURNSHIELD KITS

- Burnshield dressings provide essential physical protection and instant cooling for someone who is burnt
- Each dressing is saturated with a sterile water based gel which moisturizes the burn area, cools and soothes the skin thereby minimises trauma
- The Burnshield Responder Kit contains: 8 x BU/091, 2 x BU/092, 2 x BU/093, 3 x BU/094, 2 x BU/095, 3 x BU/098, 2 x Secure Bandages, 1 x Emergency Thermal Blanket
- The Burnshield Emergency Burncare Rescue Kit contains: 2 x BU/091, 2 x BU/092, 1 x BU/094, 1 x BU/095, 1 x BU/098, 1 x Emergency Thermal Blanket, 2 x Secure Bandages
- ISO 9001 EN 46001 CE accredited

BU/081 Burnshield Responder Kit

BU/082 Burnshield Rescue Kit



BU/081



BU/082

BURNSHIELD FIRE TRAUMA BLANKET

- A pure wool based sterile survival, Hydrogel blanket
- Provides immediate physical, heat and smoke shielding as well as relieving burn victims through cooling and trauma reducing properties
- The sterile water based gel helps smother the flames and protect the victim
- Dramatically reduces the risk of infection from 1st, 2nd or 3rd degree burns

BU/084 Fire/Trauma Blanket - 1200 x 1600mm **PRICE: £197.95**



BU/084